


# Design Primer & Probe

Course: MICROBIAL GENOME & MICROBIOME ANALYSIS - 2024

Ms. Ngo Dai Phu

PI: PhD. Luu Phuc Loi

Email: [phudaingobio@gmail.com](mailto:phudaingobio@gmail.com); [bioinformatics.mgma@gmail.com](mailto:bioinformatics.mgma@gmail.com);



### Kary B. Mullis

The Nobel Prize in Chemistry 1993

Born: 28 December 1944, Lenoir, NC, USA

Died: 7 August 2019, Newport Beach, CA, USA


Prize motivation: “for his invention of the polymerase chain reaction (PCR) method”

Prize share: ½

### Work

An organism's genome is stored inside DNA molecules, but analyzing this genetic information requires quite a large amount of DNA. In 1985, Kary Mullis invented the process known as polymerase chain reaction (PCR), in which a small amount of DNA can be copied in large quantities over a short period of time. By applying heat, the DNA molecule's two strands are separated and the DNA building blocks that have been added are bonded to each strand. With the help of the enzyme DNA polymerase, new DNA chains are formed and the process can then be repeated. PCR has been of major importance in both medical research and forensic science.

# What are primer and probe? Why do we need to design for our own?



# General Rules for Design of Oligonucleotides

With four possibilities (A, C, G, T) for selection of nucleotides at each position of the sequence,  $4^{18}$  ( $\sim 10^{11}$ ) possibilities exist for the design of an 18 nucleotide oligonucleotide.

## Lengths of PCR Primers and Products


- PCR primers in length: 18 and 24 nt
  - PCR primers of 15 nucleotides or shorter are only used for arbitrary or random short priming in the mapping of simple genomes.
- Melting temperature ( $T_m$ ): approximately  $60^\circ\text{C}$  ( $59 \pm 2^\circ\text{C}$ ) for qPCR but may vary ( $55 \pm 5^\circ\text{C}$ ) for conventional PCR
- $T_m$  difference between primers:  $\leq 5^\circ\text{C}$ 
  - To match the  $T_m$  of both forward and reverse primers, the length of a primer with high GC content needs to be shorter compared to one with less GCs.
- Primer pairs' GC content: 40–60% GC
- Product in length: approximately 200 to 1500 bp for conventional PCR and  $< 200$  bp for qPCR
- Primer dimers (hairpin, self-dimer, hetero-dimer):  $\Delta G \geq -2.0$  kcal
  - Any primer with both a terminal  $\Delta G < -2.0$  kcal and an extendable 3'-end (5'-overlap) should be avoided.
  - Avoid 3' clamping (examine the 5 bases of the 3' and accept 3 of these as A or T and 2 as G or C)
  - The strongest overall dimer should be unstable ( $\Delta G \geq -6.0$  kcal)

## Lengths of Oligonucleotide Hybridization Probes

- Hybridization probes have been applied with lengths from around 20 nt up to several hundreds of nt.
  - At least one or two nucleotide differences are needed for separation of target from nontarget sequences.
  - Mismatches near the end of probes are less destabilizing than internal mismatches
- For microarray experiments, the usual size is 25 nt, but probes of 50–70 nt are also used on certain applications

# Example 1: Design primer for Virus H5N1

- DNA sequence of the hemagglutinin (HA) gene of strain 1734 with serotype H5N1 isolated from duck in Fujian in 2005 with GenBank/EBI/DDBJ acc. no. DQ095629
- <https://www.ncbi.nlm.nih.gov/nuccore/>

 **National Library of Medicine**  
National Center for Biotechnology Information

Log in

Nucleotide

Nucleotide

Search

Advanced

Help

FASTA

### Influenza A virus (A/Duck/Fujian/1734/05(H5N1)) hemagglutinin (H5N1) cds

GenBank: DQ095629.2

[GenBank](#) [Graphics](#) [PopSet](#)

>DQ095629.2 Influenza A virus (A/Duck/Fujian/1734/05(H5N1)) hemagglutinin (HA) gene, partial cds

```
ATGGAGAAAATAGTGCTTCTTCTTGCAATAGTCAGCCTTGTTAAAAGTGATCAGATTTGCATTGGTTACC
ATGCAAACTCGACAGAGCAGGTTGACACAATAATGGAAGAAGCTTACTGTTACACATGCCCAAGA
CATACTGGAAAGACACACAACGGGAAGCTCTGCGATCTAGATGGAGTGAAGCCTCTGATTTTAAAGAGAT
TGATGTAGCTGGATGGCTCCTCGGAAACCAATGTGTGACGAATTCATCAATGTGCCGGAATGGTCTT
ACATAGTGGAGAAGGCCAACCCAGCCAATGACCTCTGTTACCCAGGGAATTTCAACGACTATGAAGAACT
GAAACACCTATTGAGCAGAATAAACCATTTTGAGAAAATTGAGATCATCCCCAAAAGTCTTGGTCCGAT
CATGAAGCCTCATCAGGGGTGAGCTCAGCATGTCCATACCAGGGAACGCCCTCCTTTTTCAGAAATGTGG
TATGGCTTATCAAAAAGAACAATACATACCAACAATAAAGAGAAGCTACAATAATACCAACGGAAGA
TCTTTTGATACTGTGGGGATTTCATCTTAATGATGCGGCAGAGCAGACAAGCTCTATCAAAACCA
ACCACCTATATTTCCGTTGGGACATCAACACTAAACCAGAGATTGGTACCAAAAATAGCTATTAGATCCA
AAGTAAACGGGCAAAGTGGAAGGATGGATTTCTTCTGGACAATTTAAAACCGAATGATGCAATCAACTT
CGAGAGTAATGGAAATTTTCATTGCTCCAGAATATGCATACAAAATTGTCAAGAAAGGGGACTCAGCAATT
ATGAAAAGTGAAGTGAATACGGTAACTGCAACACCATGTGTCAAATCCAATAGGGGCGATAAACTCTA
GTATGCCATTCCACAACATACACCTCTCACCATCGGGGAATGCCCCAAATATGTGAAATCAAACAAATT
AGTCTTGCGACTGGGCTCAGAAATAGTCCTCTAAGAGAAAGAAGAAAAAGAGGACTATTTGGAGCT
```

Send to:

☒ Complete Record

☐ Coding Sequences

☐ Gene Features

**Choose Destination**

☒ File

☐ Clipboard

☐ Collections

☐ Analysis Tool

Download 1 item.

Format

FASTA

Show GI ☐

Create File

**Related information**[Protein](#)[PubMed](#)[Taxonomy](#)

9

# Example 1: Design primer for Virus H5N1

- <https://bioinfo.ut.ee/primer3-0.4.0/>

Primer3 (v. 0.4.0) Pick primers from a DNA sequence.	<a href="#">Checks for mispriming in template.</a>	<a href="#">disclaimer</a>	<a href="#">Primer3 Home</a>
	<a href="#">Primer3plus interface</a>	<a href="#">cautions</a>	<a href="#">FAQ/WIKI</a>

There is a newer version of Primer3 available at <http://primer3.ut.ee>

Paste source sequence below (5'->3', string of ACGTNacgtn -- other letters treated as N -- numbers and blanks ignored). FASTA format ok. Please N-out undesirable sequence (vector, ALUs, LINEs, etc.) or use a [Mispriming Library \(repeat library\)](#):

```
ATGGAGAAAATAGTGCTTCTTCTTGCAATAGTCAGCCTTGTTAAAAGTGATCAGATTTGCATTGGTTACCATGCAAACAACCTCGACAGAGCAGGT
TGACACAATAATGGAAAAGAACGTTACTGTTACACATGCCCAAGACATACTGGAAAAGACACACAACGGGAAGCTCTGCGATCTAGATGGAGTGA
AGCCTCTGATTTTAAAGAGATTGTAGTGTAGCTGGATGGCTCCCGGAAACCAATGTGTGACGAATTCATCAATGTGCCGGAATGGTCTTACATA
GTGGAGAAGGCCAACCCAGCCAATGACCTCTGTACCCAGGGAATTTCAACGACTATGAAGAACTGAAACACCTATTGAGCAGAATAAACCATTT
TGAGAAAATTCAGATCATCCCCAAAAGTTCTTGGTCCGATCATGAAGCCTCATCAGGGGTGAGCTCAGCATGTCCATACCAGGGAACGCCCTCCT
TTTTCAGAAATGTGGTATGGCTTATCAAAAAGACAATACATACCCAACAATAAAGAGAAGCTACAATAATACCAACCAGGAAGATCTTTTGATA
```

<input checked="" type="checkbox"/> Pick left primer, or use left primer below:	<input type="checkbox"/> Pick hybridization probe (internal oligo), or use oligo below:	<input checked="" type="checkbox"/> Pick right primer, or use right primer below (5' to 3' on opposite strand):
<input type="text"/>	<input type="text"/>	<input type="text"/>

<a href="#">Sequence Id:</a>	<input type="text"/>	A string to identify your output.
<a href="#">Targets:</a>	<input type="text"/>	E.g. 50,2 requires primers to surround the 2 bases at positions 50 and 51. Or mark the <a href="#">source sequence</a> with [ and ]: e.g. ...ATCT[CCCC]TCAT.. means that primers must flank the central CCCC.
<a href="#">Excluded Regions:</a>	<input type="text"/>	E.g. 401,7 68,3 forbids selection of primers in the 7 bases starting at 401 and the 3 bases at 68. Or mark the <a href="#">source sequence</a> with < and >: e.g. ...ATCT<CCCC>TCAT.. forbids primers in the central CCCC.
<a href="#">Product Size Ranges</a>	<input type="text" value="100-300"/>	

<a href="#">Number To Return</a>	<input type="text" value="5"/>	<a href="#">Max 3' Stability</a>	<input type="text" value="9.0"/>
<a href="#">Max Repeat Mispriming</a>	<input type="text" value="12.00"/>	<a href="#">Pair Max Repeat Mispriming</a>	<input type="text" value="24.00"/>
<a href="#">Max Template Mispriming</a>	<input type="text" value="12.00"/>	<a href="#">Pair Max Template Mispriming</a>	<input type="text" value="24.00"/>

## General Primer Picking Conditions



# Example 1: Design primer for Virus H5N1

## Primer3 Output

No mispriming library specified  
Using 1-based sequence positions

OLIGO	start	len	tm	gc%	any	3' seq
LEFT PRIMER	697	20	60.10	45.00	2.00	0.00 TCCAAAGTAAACGGGCAAAG
RIGHT PRIMER	899	20	59.93	50.00	4.00	3.00 GCCCCTATTGGAGTTTGACA

SEQUENCE SIZE: 1698  
INCLUDED REGION SIZE: 1698

PRODUCT SIZE: 203, PAIR ANY COMPL: 5.00, PAIR 3' COMPL: 2.00

1 ATGGAGAAAATAGTGCTTCTTCTTGAATAGTCAGCCTTGTTAAAAGTGATCAGATTTGC

61 ATTGGTTACCATGCAAACAACCTCGACAGAGCAGGTTGACACAATAATGGAAAAGAACGTT

121 ACTGTTACACATGCCCAAGACATACTGGAAAAGACACACAACGGGAAGCTCTGCGATCTA

181 GATGGAGTGAAGCCTCTGATTTTAAGAGATTGTAGTGTAGCTGGATGGCTCCTCGGAAAC

241 CCAATGTGTGACGAATTCATCAATGTGCCGAATGGTCTTACATAGTGAGAAAGGCCAAC

301 CCAGCCAATGACCTCTGTTACCCAGGGAATTTCAACGACTATGAAGAACTGAAACACCTA

361 TTGAGCAGAATAAACCATTTTGAGAAAATTCAGATCATCCCCAAAAGTTCTTGGTCCGAT

421 CATGAAGCCTCATCAGGGGTGAGCTCAGCATGTCCATACCAGGGAACGCCCTCCTTTTTC

481 AGAAATGTGGTATGGCTTATCAAAAAGAACAATACATACCCAACAATAAAGAGAAGCTAC

541 AATAATACCAACCAGGAAGATCTTTTGATACTGTGGGGGATTCATCATTCTAATGATGCG

1501 ACGTATGACTACCCGAGTATTCAGAAGAAGCAAGATTA AAAAGAGAGGAAATAAGTGGA

1561 GTAAAATTGGAATCAATAGGAACCTACCAAATACTGTCAATTTATTCAACAGTGGCGAGT

1621 TCTCTAGCACTGGCAATCATGGTGGCTGGTCTATCTTTATGGATGTGCTCCAATGGGTCG

1681 TTACAATGCAGAATTGCA

KEYS (in order of precedence):

>>>>> left primer

<<<<< right primer

ADDITIONAL OLIGOS

	start	len	tm	gc%	any	3' seq
1 LEFT PRIMER	82	20	60.03	55.00	4.00	1.00 TCGACAGAGCAGGTTGACAC
RIGHT PRIMER	234	20	59.83	60.00	4.00	0.00 GAGGAGCCATCCGCTACAC
PRODUCT SIZE: 153, PAIR ANY COMPL: 5.00, PAIR 3' COMPL: 2.00						
2 LEFT PRIMER	82	20	60.03	55.00	4.00	1.00 TCGACAGAGCAGGTTGACAC
RIGHT PRIMER	282	20	59.82	50.00	4.00	2.00 GTAAGACCATTCCGGCACAT
PRODUCT SIZE: 201, PAIR ANY COMPL: 3.00, PAIR 3' COMPL: 0.00						
3 LEFT PRIMER	572	20	60.13	45.00	4.00	2.00 TGTGGGGGATTCATCATTCT
RIGHT PRIMER	716	20	60.10	45.00	2.00	0.00 CTTTGCCCGTTTACTTTGGA
PRODUCT SIZE: 145, PAIR ANY COMPL: 3.00, PAIR 3' COMPL: 1.00						
4 LEFT PRIMER	150	20	60.16	50.00	2.00	2.00 AAAGACACACAACGGGAAGC
RIGHT PRIMER	414	20	59.90	45.00	7.00	1.00 CCAAGAACTTTTGGGGATGA
PRODUCT SIZE: 265, PAIR ANY COMPL: 4.00, PAIR 3' COMPL: 0.00						

Statistics

	con	too	in	in	no	tm	tm	high	high		high	
	sid	many	tar	excl	bad	GC	too	any	3'	poly	end	
	ered	Ns	get	reg	GC%	clamp	low	high	compl	compl	X	stab
Left	14722	0	0	0	44	0	6567	3702	1	10	0	186
Right	14582	0	0	0	44	0	6424	3736	0	4	0	202

Pair Stats:

considered 790, unacceptable product size 772, ok 18

primer3 release 1.1.4

# Example 1: Design primer for Virus H5N1

- <https://www.idtdna.com/calc/analyzer/>

Left Primer – TCCAAAGTAAACGGGCAAAG  
Hairpin

idtdna.com/calc/analyzer/

Instructions | Definitions | Feedback

Sequence

5' MOD INTERNAL 3' MOD MIXED BASES

TCCAAAGTAAACGGGCAAAG

Bases 20

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default)

Target type DNA

Oligo Conc 0.25  $\mu$ M

Na<sup>+</sup> Conc 50 mM

Mg<sup>++</sup> Conc 3 mM

dNTPs Conc 0.8 mM

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER



Image Batch date: 5/29/2024 4:41 PM

TCCAAAGTAAACGGGCAAAG

Nucleotide type DNA

Na Concentration 50 mM

Mg Concentration 3 mM

Suboptimality 50 %

Sequence type Linear

Temperature 25 °C

Max Foldings 20

Start Position 0

Stop Position 0

UPDATE

ADD TO ORDER

## Homo-Dimer Analysis

The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

## Structures

structure	Image	$\Delta G$ (kcal.mole <sup>-1</sup> )	T <sub>m</sub> (°C)	$\Delta H$ (kcal.mole <sup>-1</sup> )	$\Delta S$ (cal.K <sup>-1</sup> .mole <sup>-1</sup> )	Output
1		0.24	21	-17.9	-60.86	<div>Ct</div> <div>Det</div>
2		0.45	16.3	-15.2	-52.5	<div>Ct</div> <div>Det</div>
3		1.23	0.9	-14	-51.09	<div>Ct</div> <div>Det</div>



# Example 1: Design primer for Virus H5N1

- <https://www.idtdna.com/calc/analyser/>

Left Primer – TCCAAAGTAAACGGGCAAAG  
Seft-Dimer

idtdna.com/calc/analyser/

ed Bookm... Sessions

Instructions | Definitions | Feedback

Sequence

5' MOD INTERNAL 3' MOD MIXED BASES

TCCAAAGTAAACGGGCAAAG

Bases 20

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default)

Target type DNA

Oligo Conc 0.25  $\mu$ M

Na<sup>+</sup> Conc 50 mM

Mg<sup>++</sup> Conc 3 mM

dNTPs Conc 0.8 mM

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER

## Homo-Dimer Analysis

The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

Delta G: -3.61 kcal/mol Base Pairs: 2

5' TCCAAAGTAAACGGGCAAAG

: || :

3' GAAACGGGCAAATGAAACCT

Delta G: -3.14 kcal/mol Base Pairs: 2

5' TCCAAAGTAAACGGGCAAAG

: || :

3' GAAACGGGCAAATGAAACCT

Delta G: -3.07 kcal/mol Base Pairs: 2

5' TCCAAAGTAAACGGGCAAAG

: || ::

3' GAAACGGGCAAATGAAACCT

Delta G: -3.07 kcal/mol Base Pairs: 2

5' TCCAAAGTAAACGGGCAAAG

: || :: :

3' GAAACGGGCAAATGAAACCT

Delta G: -1.34 kcal/mol Base Pairs: 2

5' TCCAAAGTAAACGGGCAAAG

: || :: :

3' GAAACGGGCAAATGAAACCT

# Example 1: Design primer for Virus H5N1

- <https://www.idtdna.com/calc/analyzer/>

Right Primer – GCCCCTATTGGAGTTTGACA  
Hairpin

## OligoAnalyzer

Sequence

5' MOD ▾ INTERNAL ▾ 3' MOD ▾ MIXED BASES ▾

GCCCCTATTGGAGTTTGACA

Bases 20

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default) ▾

Target type

DNA ▾

Oligo Conc

0.25

μM

Na<sup>+</sup> Conc

50

mM

Mg<sup>++</sup> Conc

3

mM

dNTPs Conc

0.8

mM

Instructions | Definitions | Feedback

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER

Image Batch date: 5/29/2024 4:37 PM

GCCCCTATTGGAGTTTGACA

Nucleotide type

DNA ▾

Sequence type

Linear ▾

Na Concentration

50

mM

Temperature

25

°C

Mg Concentration

3

mM

Max Foldings

20

Suboptimality

50

%

Start Position

0

Stop Position

0

UPDATE

ADD TO ORDER

## Hetero-Dimer Analysis

The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

## Structures

structure	Image	ΔG (kcal.mole <sup>-1</sup> )	T <sub>m</sub> (°C)	ΔH (kcal.mole <sup>-1</sup> )	ΔS (cal.K <sup>-1</sup> .mole <sup>-1</sup> )	Output
1		-1.1	42.5	-19.8	-62.73	<div>Ct</div> <div>Det</div>
2		-0.99	43.5	-16.9	-53.38	<div>Ct</div> <div>Det</div>
3		-0.53	34.3	-17.6	-57.24	<div>Ct</div> <div>Det</div>

10

# Example 1: Design primer for Virus H5N1

- <https://www.idtdna.com/calc/analyser/>

Right Primer – GCCCCTATTGGAGTTTGACA  
Seft-Dimer

## OligoAnalyzer

Sequence

5' MOD INTERNAL 3' MOD MIXED BASES

GCCCCTATTGGAGTTTGACA

Bases 20

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default)

Target type DNA

Oligo Conc 0.25 μM

Na<sup>+</sup> Conc 50 mM

Mg<sup>++</sup> Conc 3 mM

dNTPs Conc 0.8 mM

Instructions | Definitions | Feedback

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER



## Hetero-Dimer Analysis

The delta G is calculated by taking into account the longest stretch of complementary bases. These pairs of complementary bases are represented by a solid line. Dotted lines represent additional complementary bases for that dimer structure, but their presence does not impact calculated delta G values. Actual delta G values may vary based on presence of additional complementary bases. The Maximum Delta G value refers to the free energy of the oligo sequence binding to its perfect complement.

Delta G: -3.14 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -1.95 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -3.07 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -1.95 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -3.07 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -1.6 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -3.07 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

Delta G: -1.34 kcal/mol Base Pairs: 2

5' GCCCCTATTGGAGTTTGACA

3' ACAGTTTGAGGTTATCCCCG

# Example 1: Design primer for Virus H5N1

- <https://www.idtdna.com/calc/analyzer/>

Left Primer – TCCAAAGTAAACGGGCAAAG

Right Primer – GCCCCTATTGGAGTTTGACA

## Hetero-Dimer

Sequence

5' MOD ▾ INTERNAL ▾ 3' MOD ▾ MIXED BASES ▾

TCC AAA GTA AAC GGG CAA AG

Bases 20

CLEAR SEQUENCE

Try the new batch mode here

Parameter sets

SpecSheet (Default) ▾

Target type

DNA ▾

Oligo Conc

0.25

μM

Na<sup>+</sup> Conc

50

mM

Mg<sup>2+</sup> Conc

3

mM

dNTPs Conc

0.8

mM

ANALYZE

HAIRPIN

SELF-DIMER

HETERO-DIMER

NCBI BLAST

TM MISMATCH

ADD TO ORDER

## Hetero-Dimer Analysis

Primary Sequence:

5'- TCC AAA GTA AAC GGG CAA AG

Secondary Sequence:

5' GCCCCTATTGGAGTTTGACA 3'

CREATE COMPLEMENT CALCULATE



# General Rules for Design of Oligonucleotides

With four possibilities (A, C, G, T) for selection of nucleotides at each position of the sequence,  $4^{18}$  ( $\sim 10^{11}$ ) possibilities exist for the design of an 18 nucleotide oligonucleotide.

## Lengths of PCR Primers and Products

- PCR primers in length: 18 and 24 nt
  - PCR primers of 15 nucleotides or shorter are only used for arbitrary or random short priming in the mapping of simple genomes.
- Melting temperature ( $T_m$ ): approximately  $60\text{ }^{\circ}\text{C}$  ( $59\pm 2\text{ }^{\circ}\text{C}$ ) for qPCR but may vary ( $55\pm 5\text{ }^{\circ}\text{C}$ ) for conventional PCR
- $T_m$  difference between primers:  $\leq 5\text{ }^{\circ}\text{C}$ 
  - To match the  $T_m$  of both forward and reverse primers, the length of a primer with high GC content needs to be shorter compared to one with less GCs.
- Primer pairs' GC content: 40–60% GC
- Product in length: approximately 200 to 1500 bp for conventional PCR and  $< 200$  bp for qPCR
- Primer dimers (hairpin, self-dimer, hetero-dimer):  $\Delta G \geq -2.0\text{ kcal}$ 
  - Any primer with both a terminal  $\Delta G < -2.0\text{ kcal}$  and an extendable 3'-end (5'-overlap) should be avoided.
  - Avoid 3' clamping (examine the 5 bases of the 3' and accept 3 of these as A or T and 2 as G or C)
  - The strongest overall dimer should be unstable ( $\Delta G \geq -6.0\text{ kcal}$ )

## Lengths of Oligonucleotide Hybridization Probes

- Hybridization probes have been applied with lengths from around 20 nt up to several hundreds of nt.
  - At least one or two nucleotide differences are needed for separation of target from nontarget sequences.
  - Mismatches near the end of probes are less destabilizing than internal mismatches
- For microarray experiments, the usual size is 25 nt, but probes of 50–70 nt are also used on certain applications


# Homework 1

1. Between global and local alignments, which algorithm do we use for designing primer/probe?
2. Is the primer pair appropriate for the hemagglutinin (HA) gene of strain 1734 with serotype H5N1 and Why?

Based on '**General Rules for Design of Oligonucleotides**' to answer the question 2.

# Example 2: Design primer for Virus H7N2

The serotype H7N2 of influenza virus: acc. no. U20461

 **National Library of Medicine**  
National Center for Biotechnology Information

Log in

Nucleotide Nucleotide  Search

Advanced Help

GenBank Send to: Change region shown

## Influenza A virus (A/duck/Hong Kong/293/78(H7N2)) hemagglutinin mRNA, complete cds

GenBank: U20461.1

[FASTA](#) [Graphics](#)

[Go to:](#) ☒

LOCUS	IAU20461	1731 bp	mRNA	linear	VRL 10-MAR-2006
DEFINITION	Influenza A virus (A/duck/Hong Kong/293/78(H7N2)) hemagglutinin precursor (HA) mRNA, complete cds.				
ACCESSION	U20461				
VERSION	U20461.1				
KEYWORDS	.				
SOURCE	Influenza A virus (A/duck/Hong Kong/293/1978(H7N2))				
ORGANISM	<a href="#">Influenza A virus (A/duck/Hong Kong/293/1978(H7N2)).</a> Viruses; Riboviria; Orthornavirae; Negarnaviricota; Polyploviricota; Insthoviricetes; Articulavirales; Orthomyxoviridae; Alphainfluenzavirus; Alphainfluenzavirus influenzae.				
REFERENCE	1 (bases 1 to 1731)				
AUTHORS	Rohm,C., Horimoto,T., Kawaoka,Y., Suss,J. and Webster,R.G.				
TITLE	Do hemagglutinin genes of highly pathogenic avian influenza viruses constitute unique phylogenetic lineages?				
JOURNAL	Virology 209 (2), 664-670 (1995)				
PUBMED	<a href="#">7778300</a>				
REFERENCE	2 (bases 1 to 1731)				
AUTHORS	Rohm,C.				

☒ Complete Record  
☐ Coding Sequences  
☐ Gene Features

**Choose Destination**  
☒ File ☐ Clipboard  
☐ Collections ☐ Analysis Tool

Download 1 item.  
Format  
FASTA  
Show GI ☐  
Create File

Retrieve, view, and download influenza virus genomic and protein sequences.

**Related information**

Protein

PubMed

Taxonomy

Full text in PMC

15



<https://www.ncbi.nlm.nih.gov/tools/primer-blast/>



## Primer-BLAST

### Finding primers specific to your PCR template (using Primer3 and BLAST).

Primers common for a group of sequences

[Retrieve recent results](#)   [Publication](#)   [Tips for finding specific primers](#)

Save search parameters

[Reset page](#)

Clear

Range ?

Clear

	From	To
Forward primer		
Reverse primer		

Or, upload FASTA file

Choose File No file chosen

Use my own forward primer  
(5'→3' on plus strand)

?

Use my own reverse primer (5'-  
→3' on minus strand)

?

PCR product size

Min	Max
70	1000

# of primers to return

10

### Primer melting temperatures ( $T_m$ )

Min	Opt	Max	Max $T_m$ difference
57.0	60.0	63.0	3

# Example 2: Design primer for Virus H7N2

Note: Parameter values that differ from the default are highlighted in yellow

## Primer Pair Specificity Checking Parameters

### Specificity check

☒ Enable search for primer pairs specific to the intended PCR template ?

### Search mode

User guided ?

### Database

nr ?

### Exclusion

☐ Exclude predicted Refseq transcripts (accession with XM, XR prefix) ☐ Exclude uncultured/environmental sample sequences ?

### Organism

Influenzavirus A

Add organism

Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type. ?

### Entrez query (optional)

?

### Primer specificity stringency

Primer must have at least 2 total mismatches to unintended targets, including at least 2 mismatches within the last 5 bps at the 3' end. ?

Ignore targets that have 6 or more mismatches to the primer. ?

### Max target amplicon size

4000 ?

### Allow splice variants

☐ Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input) ?

Get Primers

☒ Show results in a new window ☒ Use new graphic view ?

# Example 2: Design primer for Virus H7N2

## Primer-BLAST

A tool for finding specific primers

Finding primers specific to your PCR template (using Primer3 and BLAST)

Input PCR template lcl|Query\_1

Range 1 - 1731

Your PCR template is highly similar to the following sequence(s) from the search database. To increase the chance of finding specific primers, please review the list below and select all sequences (within the given sequence ranges) that are intended or allowed targets.

Select: All None Selected:9

Accession	Title	Identity	Alignment length	Seq. start	Seq. stop
<input type="checkbox"/> U20461.1	Influenza A virus (A/duck/Hong Kong/293/78(H7N2)) hemagglutinin precursor (HA) mRNA, complete cds	100%	1731	1	1731
<input checked="" type="checkbox"/> CY006029.1	Influenza A virus (A/dk/Hong Kong/293/1978(H7N2)) segment 4, complete sequence	99.94%	1731	2	1732
<input checked="" type="checkbox"/> CY067686.1	Influenza A virus (A/swine/KU/16/2001(H7N2)) segment 4 sequence	99.36%	1731	2	1732
<input checked="" type="checkbox"/> AB302789.1	Influenza A virus (A/duck/Hong Kong/301/1978(H7N2)) HA gene for haemagglutinin, complete cds	99.82%	1698	1	1698
<input checked="" type="checkbox"/> KX130802.1	Influenza A virus (A/twite/Gangcha/02/2006(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.76%	1683	1	1683
<input checked="" type="checkbox"/> KX130803.1	Influenza A virus (A/tree sparrow/Gangcha/02/2006(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.7%	1683	1	1683
<input checked="" type="checkbox"/> KX130801.1	Influenza A virus (A/twite/Gangcha/01/2006(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.7%	1683	1	1683
<input checked="" type="checkbox"/> KX130820.1	Influenza A virus (A/pika/Maduo/01/2009(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.64%	1683	1	1683
<input checked="" type="checkbox"/> KX130811.1	Influenza A virus (A/pika/QH-Maduo/02/2006(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.64%	1683	1	1683
<input checked="" type="checkbox"/> KX130809.1	Influenza A virus (A/ground jay/Maduo/01/2006(H7N2)) segment 4 hemagglutinin (HA) gene, complete cds	99.64%	1683	1	1683
<input type="checkbox"/> DQ003216.1	Influenza A virus (A/duck/Hongkong/301/72(H7N1)) hemagglutinin (HA) gene, complete cds	99.58%	1682	1	1682

Submit

☐

Show results in a new window


# Example 2: Design primer for Virus H7N2

Input PCR template Id|Query\_1  
Range 1 - 1731

Specificity of primers Primers may **not** be specific to the input PCR template as targets were found in selected database:Nucleotide collection (nt) (Organism limited to Alphainfluenzavirus)...[help on specific primers](#)

Other reports [▶Search Summary](#)

## — Graphical view of primer pairs



## — Detailed primer reports

You can re-search for specific primers by accepting some of the unintended targets, check the box(es) next to the ones you accept and try again to re-search for specific primers [Submit](#) [?](#)

### Primer pair 1

	Sequence (5'->3')	Template strand	Length	Start	Stop	Tm	GC%	Self complementarity	Self 3' complementarity
Forward primer	GTCACCTGGACCACCTCAGTG	Plus	20	252	271	59.97	60.00	5.00	3.00
Reverse primer	ATCGATCCGTCAGATTGGC	Minus	20	740	721	59.97	55.00	6.00	2.00
Product length	489								

## Homework 2

1. Are all primers of 'Products of intended targets' for Primer pair 1 matching the type H7N2?
2. Are all 'Products of unintended targets' forming mismatches for Primer pair 1 that are expected not to result in generation of product?

Based on 'General Rules for Design of Oligonucleotides' to answer this question.

# Reference

- <https://www.sigmaaldrich.com/NL/en/technical-documents/technical-article/genomics/pcr/assay-optimization-and-validation>
- <https://www.sigmaaldrich.com/NL/en/technical-documents/protocol/genomics/pcr/pcr-qpcr-dpcr-assay-design>
- <https://www.idtdna.com/pages/education/decoded/article/designing-pcr-primers-and-probes>
- <https://www.idtdna.com/pages/education/decoded/article/design-efficient-pcr-and-qpcr-primers-and-probes-using-online-tools>
- Christensen, Henrik, and John Elmerdahl Olsen. "Primer Design: Design of Oligonucleotide PCR Primers and Hybridization Probes." *Introduction to Bioinformatics in Microbiology* (2018): 81-102.
- Nolan, Tania, Jim Huggett, and Elena Sanchez. "Good practice guide for the application of quantitative PCR (qPCR)." Teddington: LGC (2013).


# T<sub>m</sub> Calculations

- Estimation of T<sub>m</sub> by Formula
  - For oligonucleotides of ≤ 20 nt:  $T_m = 4 (G + C + 2(A-T))$
  - For oligonucleotides of 10–50 nt:  $T_m = 81.5 + 16.6\log M + 0.41(G + C\%) - (820/(\text{length of probe}))$
  - For oligonucleotides of ≥ 50 nt:  $T_m = 81.5^\circ\text{C} + 16.6\log M + 0.41(G + C\%) - (500/\text{length of probe})$
- Formamide Considerations
  - 1% formamide reduces T<sub>m</sub> by 0.6°C
- Estimation of T<sub>m</sub> by Nearest Neighbor Prediction
  - $T_d = \Delta H/(\Delta S + R\ln C) + 16.6\log_{10} [M]$
- Combination of Formula and nearest neighbor calculations
  - $T_a = 0.3 T_m \text{ primer} + 0.7 T_m \text{ product} - 14.9$




# Design of Primers for PCR and “Kwok’s Rules”

Prediction of amplication from interactions between PCR oligonucleotide primer at the 3 end and template (based on Kwok et al. 1990). Low degree of mismatch is assumed between oligonucleotide and template for the remaining part of sequence



# Design of Probes for Hybridization

Prediction of binding in hybridization between probe and template. A low degree of mismatch formation is assumed between oligonucleotide and template for the remaining part of the sequence



Program	Description	URL and reference
Amplifx2.1	PCR design program for MacOS X	<a href="https://macdownload.informer.com/amplifx/">https://macdownload.informer.com/amplifx/</a>
FastPCR	PCR primer design including multiplex PCR	<a href="http://primerdigital.com/fastpcr.html">http://primerdigital.com/fastpcr.html</a>
Netprimer	Numerous parameters of single PCR primers	<a href="http://www.premierbiosoft.com/netprimer/index.html">http://www.premierbiosoft.com/netprimer/index.html</a>
Oligo	Nearest neighbor calculations of secondary structures and T <sub>m</sub>	<a href="https://www.oligo.net/">https://www.oligo.net/</a> Rychlik and Rhoads (1989)
OligoAnalyzer	Nearest neighbor parameters. Hair-pin and primer-dimer analysis	<a href="https://eu.idtdna.com/pages/tools">https://eu.idtdna.com/pages/tools</a> Owczarzy et al. (2008)
OligoCalc	Physical properties of oligonucleotides, self complementarity and hairpin loop formation	<a href="http://biotools.nubic.northwestern.edu/OligoCalc.html">http://biotools.nubic.northwestern.edu/OligoCalc.html</a> Kibbe (2007)
OligoFaktory	DNA microarrays, primers for PCR, siRNAs	<a href="http://www.bioinformatics.org/oligofaktory/">http://www.bioinformatics.org/oligofaktory/</a>
Pride and Genome pride	PCR and microarray	<a href="http://pride.molgen.mpg.de/">http://pride.molgen.mpg.de/</a> (Staden package)
Pria	Search for primers in multiple alignments	<a href="https://services.birc.au.dk/pria/">https://services.birc.au.dk/pria/</a> Fredslund et al. (2005)
Primegens	Primer design	<a href="http://primegens.org/">http://primegens.org/</a> Xu et al. (2002)
Primer3	The most frequently used program see also Activity 5.1	<a href="http://bioinfo.ut.ee/primer3-0.4.0/primer3/input.htm">http://bioinfo.ut.ee/primer3-0.4.0/primer3/input.htm</a> Rozen and Skaletsky (2000)
Primer3Plus	Simplified version of Primer3	<a href="http://www.bioinformatics.nl/cgi-bin/primer3plus/primer3plus.cgi/">http://www.bioinformatics.nl/cgi-bin/primer3plus/primer3plus.cgi/</a>
Primer Design Assistant (PDA)	Primer design with all parameters	<a href="http://dbb.nhri.org.tw/primer/">http://dbb.nhri.org.tw/primer/</a> Chen et al. (2003)
PrimerPremier	PCR primer design, multiplex, degeneracies and nested PCRs	<a href="http://www.premierbiosoft.com/primerdesign/index.html">http://www.premierbiosoft.com/primerdesign/index.html</a>
PrimerQuest	PCR and hybridization probes. Based on Primer3	<a href="https://eu.idtdna.com/PrimerQuest/Home/Index">https://eu.idtdna.com/PrimerQuest/Home/Index</a>
Primer Search (EMBOSS)	Matches between primers-pairs and DNA template by string comparisons	<a href="http://emboss.sourceforge.net/">http://emboss.sourceforge.net/</a> Rice et al. (2000)
Primo Oligo	Calculation of T <sub>m</sub>	<a href="http://www.changbioscience.com/primo/oligo.html">http://www.changbioscience.com/primo/oligo.html</a>
Web Primer	PCR primer design	<a href="https://www.yeastgenome.org/cgi-bin/web-primer">https://www.yeastgenome.org/cgi-bin/web-primer</a>
SciTools	A series of tools	<a href="http://www.idtdna.com/SciTools/Scitools.aspx">http://www.idtdna.com/SciTools/Scitools.aspx</a> Owczarzy et al. (2008)

Program	Description	URL and reference
ARB	Design of rRNA targeted probes	<a href="http://www.arb-home.de">http://www.arb-home.de</a> Ludwig et al. (2004).
ArrayDesigner	Oligo-arrays, cDNA arrays and SNP arrays	<a href="http://www.premierbiosoft.com/dnamicroarray/index.html">http://www.premierbiosoft.com/dnamicroarray/index.html</a>
Assembly PCR oligo maker	PCR-based construction of long DNA molecules for RNA molecules by T7 RNA polymerase	<a href="http://www.yorku.ca/pjohnson/AssemblyPCRoligomaker.html">http://www.yorku.ca/pjohnson/AssemblyPCRoligomaker.html</a>
Beacon designer™	Real-time PCR	<a href="http://www.premierbiosoft.com/molecular_beacons/index.html">http://www.premierbiosoft.com/molecular_beacons/index.html</a>
BLOCKMAKER and CODEHOP	Degenerate primers to genes of proteins	<a href="https://4virology.net/virology-ca-tools/j-codehop/">https://4virology.net/virology-ca-tools/j-codehop/</a> Rose et al. (2003)
Expeditior	QTL design	<a href="https://www.animalgenome.org/cgi-bin/expeditor/expeditor2">https://www.animalgenome.org/cgi-bin/expeditor/expeditor2</a>
Geneasher2	Degenerate PCR primers based on multiple aligned sequences	<a href="https://bibiserv.cebitec.uni-bielefeld.de/geneasher2/">https://bibiserv.cebitec.uni-bielefeld.de/geneasher2/</a>
MFOLD	Evaluates probes	<a href="http://www.unafold.org/mfold/applications/rna-folding-form-v2.php">http://www.unafold.org/mfold/applications/rna-folding-form-v2.php</a>
OligoPicker	microarray design of oligonucleotides	<a href="https://pga.mgh.harvard.edu/oligopicker/index.html">https://pga.mgh.harvard.edu/oligopicker/index.html</a>
Pira-PCR	SNP	<a href="http://primer1.soton.ac.uk/primer2.html">http://primer1.soton.ac.uk/primer2.html</a> Ke et al. (2001)
PrimerD	Degenerate primer pairs	<a href="http://mblab.wustl.edu/software.html#primerdLink">http://mblab.wustl.edu/software.html#primerdLink</a>
Primer exmplore	LAMP primers	<a href="http://primerexplorer.jp/e/">http://primerexplorer.jp/e/</a>
PrimerX	Site-directed mutagenesis	<a href="http://www.bioinformatics.org/primerx/cgi-bin/DNA_1.cgi">http://www.bioinformatics.org/primerx/cgi-bin/DNA_1.cgi</a>
ProDesign	Oligonucleotide design for microarray	<a href="http://wwwlabs.uhnresearch.ca/tillier/ProDesign/ProDesign.html">http://wwwlabs.uhnresearch.ca/tillier/ProDesign/ProDesign.html</a> Feng and Tillier (2007)